

## Catalog Number

48310-20

**Tubing, PFA, 0.020 inch (0.50 mm) ID, 1/16th inch (1.6 mm) OD, low pressure, 5 meter roll, ARE-Applied Research brand. 1 EA.**



### Additional Info

This PFA (*Perfluoroalkyl*) low pressure tubing is very similar to Teflon but is less permeable, more transparent and performs better in higher temperatures and pressures.

Do not use this tubing with heat over 260°C.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Color: Natural, Translucent
- Wall Thickness: 0.55 mm

**Price**  
**\$33.70**

## Additional Info

- Pressure rating: 1450 psi (99.9 bar)

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$33.70**