


Catalog Number

48410-10-C

Tubing, FEP, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD, low pressure, 25 meter roll, ARE-Applied Research brand 1 EA.



Additional Info

This FEP (*Fluoroethylene Propylene*) low pressure tubing is chosen for use in Ion Chromatography Applications, but can be used in all low-pressure, laboratory applications including before and after an HPLC Pump such as with Mobile Phase Delivery and Solvent Waste Systems. Also known as Teflon FEP.

Do not use this tubing with heat over 80°C.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Color: Natural, Translucent
- Wall Thickness: 0.675 mm
- Pressure rating: 1650 psi (113.8 bar)
- Supplied on a convenient spool.

Other lengths available:

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$361.13

Printed On 2024/05/18