


# Catalog Number

48445-65-1M

**Tubing, FEP, 1/16th inch (1.6 mm) ID, 1/8th inch (3.2 mm) OD, low pressure, 1 meter roll, ARE-Applied Research brand. 1 EA.**



## Additional Info

This FEP (*Fluoroethylene Propylene*) low pressure tubing is chosen for use in Ion Chromatography Applications, but can be used in all low-pressure, laboratory applications including before and after an HPLC Pump such as with Mobile Phase Delivery and Solvent Waste Systems. Also known as Teflon FEP.

Do not use this tubing with heat over 80°C.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Color: Natural, Translucent

**Price**  
**\$18.88**

## Additional Info

- Wall Thickness: 0.80 mm
- Pressure rating: 1000 psi (68.9 bar)

NOTE: *This item is in compliance with the RoHS EU Directive 2011/65/EU and its amendment (EU) 2015/863 regarding the restriction of hazardous substances.*

Other lengths available:

- [3 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

## MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$18.88**