


## Catalog Number

48445-65-3M

**Tubing, FEP, 1/16th inch (1.6 mm) ID, 1/8th inch (3.2 mm) OD, low pressure, 3 meter roll, ARE-Applied Research brand. 1 EA.**



### Additional Info

This FEP (*Fluoroethylene Propylene*) low pressure tubing is chosen for use in Ion Chromatography Applications, but can be used in all low-pressure, laboratory applications including before and after an HPLC Pump such as with Mobile Phase Delivery and Solvent Waste Systems. Also known as Teflon FEP.

Do not use this tubing with heat over 80°C.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Color: Natural, Translucent
- Wall Thickness: 0.80 mm
- Pressure rating: 1000 psi (68.9 bar)

Other lengths available:

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$43.34**

Printed On 2024/05/06