

Catalog Number

48510-30-10M

Tubing, ETFE, 0.030 inch (0.75 mm) ID, 1/16th inch (1.6 mm) OD, low pressure, 10 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This ETFE (*Tefzel*) tubing, is extruded from virgin ETFE without any manufacturing additives. This tubing is used widely in all low-pressure, laboratory applications including before and after an HPLC Pump such as with Mobile Phase Delivery and Solvent Waste Systems. ETFE tubing is chemically inert except for some halogenated compounds and has low permeability compared with PTFE tubing.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Color: Natural, Translucent
- Wall Thickness: 0.425 mm
- Pressure rating: 1625 psi (bar 112.1)

Other lengths available:

- [1 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$70.20

Printed On 2024/04/28