


## Catalog Number

48555-65-C

**Tubing, ETFE, 1/16th inch (1.6 mm) ID, 1/8th inch (3.2 mm) OD, low pressure, 25 meter roll, ARE-Applied Research brand. 1 EA.**



### Additional Info

This ETFE (*Tefzel*) tubing, is extruded from virgin ETFE without any manufacturing additives. This tubing is used widely in all low-pressure, laboratory applications including before and after an HPLC Pump such as with Mobile Phase Delivery and Solvent Waste Systems. ETFE tubing is chemically inert except for some halogenated compounds and has low permeability compared with PTFE tubing.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Color: Natural, Translucent
- Wall Thickness: 0.80 mm
- Pressure rating: 1000 psi (bar 68.9)
- Supplied on a convenient spool.

Other lengths available:

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$385.20**

Printed On 2024/05/05