

## Catalog Number

49025-01N

**Plug. PEEK, 1/4-28, ARE-Applied Research brand. 1 EA.**



### Additional Info

Knurled for a tight grip and also has "flat" surfaces for using a wrench to quarter tighten for best seal or to loosen if plug gets stuck.

Screw thread for 1/4 - 28" holes.

Use with Delivery Cap [58820-06](#), to prevent oxygen from entering the bottle while one of the holes are not in use.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$8.48**

Printed On 2024/02/25