

## Catalog Number

49027-03Y

**End Plugs**, for HPLC Columns. Nylon 10-32 Thread. Yellow. ARE-Applied Research brand. 10/EA.



### Additional Info

Use these nylon, endstopper fittings to close and seal any guard, analytical or prep HPLC column with 10-32 threads and prevent it from drying out or from solvent leaking from the column. *Compatible with any brand of column.*

Use when storing columns or when columns are not attached to instruments will increase the column lifetime.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$8.66**

Printed On 2024/04/28