

Catalog Number

49027-03Y

End Plugs, for HPLC Columns. Nylon 10-32 Thread. Yellow. ARE-
Applied Research brand. 10/PK.



Additional Info

Use these nylon, endstopper fittings to close and seal any guard, analytical or prep HPLC column with 10-32 threads and prevent it from drying out or from solvent leaking from the column. *Compatible with any brand of column.*

Use when storing columns or when columns are not attached to instruments will increase the column lifetime.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$8.66

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$8.66

Printed On 2024/08/06