

Catalog Number

49029-02-010

**Check Valve**, luer, pocket for 0.122-inch (3.1 mm) OD tubing to male locking luer. Clear-Transparent MABS and Blue-Transparent SAN with Silicone Diaphragm. 10/PK.



### Additional Info

**This** Check Valve is a high-performance, luer valve designed for secure and reliable fluid control in laboratory applications. This valve features a clear-transparent MABS body and a blue-transparent SAN cap, with a silicone diaphragm for optimal performance.

### Key Features:

- **Material:** Constructed from medical-grade materials, including clear-transparent MABS (Methyl Methacrylate Acrylonitrile Butadiene Styrene) and blue-

**Price**  
**\$32.41**

## Additional Info

transparent SAN (Styrene Acrylonitrile), ensuring durability and chemical resistance.

- **Silicone Diaphragm:** Equipped with a high-quality silicone diaphragm that provides a reliable seal and prevents backflow.
- **Design:** Features a pocket for 0.122-inch (3.1 mm) OD tubing, allowing for easy and secure connections to male locking luer fittings.
- **Visibility:** The clear and blue-transparent materials enhance visibility, making it easier to inspect fluid flow and connections.

## Benefits:

- **Leak-Proof:** The precise design and high-quality materials ensure a secure, leak-proof connection, critical for maintaining fluid integrity in medical and laboratory settings.
- **Versatile:** Suitable for a wide range of applications, including infusion, drainage, and irrigation systems.
- **Reliable:** Manufactured to stringent quality standards, ensuring consistent performance and reliability in critical applications.

## Applications:

- **Laboratory Instruments:** Suitable for use in analytical instruments, fluid transfer systems, and other laboratory equipment where precise fluid control is essential.

## MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Price

**\$32.41**