

Catalog Number

49101-10

PQ Kit, Complete. 150 mm HPLC column, test solutions, accessories & software for performance qualification. Chemical Solutions brand. 1 EA.

**Additional Info**

This complete PQ Kit™ for self-compliance of HPLC performance qualification (PQ) provides everything you need to complete a full qualification of your HPLC instruments and perform holistic and modular AIQ tests. You can do preventative maintenance with immediate instrument uptime without waiting for vendors to verify your instrument. Fulfill your requirements to “maintain instruments in a qualified state.” See USP <1058>. You can also use this kit to qualify your HPLC instruments for accuracy and precision.

This PQ Kit includes:

Price**\$2,699.72**

Additional Info

- 6 Linearity Solutions - 30 ml (L2 and L3 60 ml)
- 1 Resolution Test Mixture - 30 ml
- 1 Wavelength Calibration Solution (Holmium Oxide) - 4 ml
- 1 Gradient Visualization Solution - 30 ml
- 1 Cogent HPLC Certified Test Column - 150 mm x 4.6 mm
- 1 Software Template. Excel with Operator Determined Criteria
- 1 Accessories for Flow Cell Filling Pump Kit
- 1 Quick Start Instructions

Click [HERE](#) for specifications of solutions, Click [HERE](#) for more information about the HSQ kit. With proper use, you can perform approximately 25 qualifications per kit.

Solutions will be supplied with at least six months until the expiration date. They are designed to be used promptly and often. Expiration date is two years after the date of manufacture.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,699.72