

## Catalog Number

49101-25S

**Performance Qualification Solutions.** Replacement set of solutions for use with the HSQ or PQ kits. Chemical Solutions brand. 1/EA



### Additional Info

This set of NIST traceable solutions are replacements for use with the HSQ™ or PQ™ kits for performance qualification of HPLC instruments for accuracy and precision.

Set includes:

- 6 Linearity Solutions - 30 ml (L2 and L3 60 ml)
- 1 Resolution Test Mixture - 30 ml
- 1 Wavelength Calibration Solution (Holmium Oxide) - 4 ml
- 1 Gradient Visualization Solution - 30 ml

Click [HERE](#) for specifications and concentrations.

*Solutions will be supplied with at least six months until the expiration date. They are designed to be used promptly and often. Expiration date is two years after the date of manufacture.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,329.96**

Printed On 2024/05/07