

## Catalog Number

49115-25

**UV-vis Calibration Solutions.** Replacement set of solutions for use with UV-vis calibration kit. Chemical Solutions brand. 1 EA.



### Additional Info

This set of NIST traceable solutions are replacements for use with the Chemical Solutions UV-vis Calibration Kit™ for performance qualification (PQ) of UV or UV-vis spectrophotometers and detectors.

Set includes:

- Absorbance Solutions - 7 x 30 ml each
- Reference Solution Diluent - 60 ml
- Wavelength Calibration Solution - 30 ml

Price

**\$1,269.42**

## Additional Info

- Stray Light Solution - 30 ml
- Certificate of Analysis.

Click [HERE](#) for more information about the solutions and the kit or Click [HERE](#) for specifications and concentrations.

*Solutions will be supplied with at least six months until the expiration date. They are designed to be used promptly and often. Expiration date is two years after the date of manufacture.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,269.42**