



Catalog Number

49305-07-3M

Tubing, PEEK, 0.007 inch (0.18 mm) ID, 1/16th inch (1.6 mm) OD, general grade, solid yellow, 3 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This general grade tubing is extruded from virgin PEEK material without any manufacturing additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.71 mm
- Color code: yellow
- Pressure rating: 5800 psi (399.9 bar)

Other lengths available:

Price
\$61.55



Additional Info

- [1 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)

Custom cuts also available. Please contact [Customer Service](#).

NOTE: *This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.*

MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$61.55