


Catalog Number

49305-20

Tubing, PEEK, 0.020 inch (0.50 mm) ID, 1/16th inch (1.6 mm) OD, general grade, solid orange, 5 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This general grade tubing is extruded from virgin PEEK material without any manufacturing additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.55 mm
- Color code: orange
- Pressure rating: 4500 psi (310.3 bar)

Other lengths available:

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$96.30

Printed On 2024/08/06

Additional Info

- [1 meter roll](#)
- [3 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

NOTE: This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$96.30

Printed On 2024/08/06