


## Catalog Number

49305-20-3M

**Tubing, PEEK, 0.020 inch (0.50 mm) ID, 1/16th inch (1.6 mm) OD, general grade, solid orange, 3 meter roll, ARE-Applied Research brand. 1 EA.**



### Additional Info

This general grade tubing is extruded from virgin PEEK material without any manufacturing additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.55 mm
- Color code: orange
- Pressure rating: 4500 psi (310.3 bar)

Other lengths available:

**Price**  
**\$59.08**

## Additional Info

- [1 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)

Custom cuts also available. Please contact [Customer Service](#).

**NOTE:** *This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$59.08**