


Catalog Number

49305-504G

Tubing, PEEK, Coiled, 0.030 (0.75 mm) ID, 1/16th inch (1.6 mm) OD.
710 mm length, HPLC grade, green stripe, with tubing connectors. ARE-
Applied Research brand. 1 EA.



Additional Info

These Flexible, Coiled PEEK Tubes will expand and contract which will save space in places such as behind an instrument and provide a Pulse Dampener Effect for Parabolic Flow Fluids which will result in a quieter base line and also reduce band broadening due to excess tubing lengths.

This coil is 10 mm wide when resting and 30 mm when fully extended, with a total tubing length of 710 mm, end to end.

Type 2 FlexiTube, color coded: green stripe for 0.030 inch (0.75 mm) ID

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$55.85

Printed On 2024/05/05