


Catalog Number

49306-501SB

Tubing, PEEK, Coiled, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD. 945 mm L, HPLC grade, blue stripe, with tubing connectors. ARE-Applied Research brand. 1 EA.



Additional Info

These Flexible, Coiled PEEK Tubes will expand and contract which will save space in places such as behind an instrument and provide a Pulse Dampener Effect for Parabolic Flow Fluids which will result in a quieter base line and also reduce band broadening due to excess tubing lengths.

This coil is 20 mm wide when resting and 120 mm when fully extended, with a total tubing length of 945 mm, end to end.

Type 3 FlexiTube, color coded: stripe blue for 0.010 inch (0.25 mm) ID

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$57.78

Printed On 2024/05/14