


# Catalog Number

49308-503Y

**Tubing, PEEK, Coiled, 0.007 (0.18 mm) ID, 1/16th inch (1.6 mm) OD.**  
1585 mm length, HPLC grade, solid yellow, with tubing connectors. ARE-  
Applied Research brand. 1 EA.



## Additional Info

These Flexible, Coiled PEEK Tubes will expand and contract which will save space in places such as behind an instrument and provide a Pulse Dampener Effect for Parabolic Flow Fluids which will result in a quieter base line and also reduce band broadening due to excess tubing lengths.

This coil is 50 mm wide when resting and 300 mm when fully extended, with a total tubing length of 1585 mm, end to end.

*Type 5 FlexiTube, color coded: yellow for 0.007 inch (0.18 mm) ID*

**Price**  
**\$63.56**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$63.56**