


Catalog Number

49310-501B

Tubing, PEEK, Coiled, 0.010 (0.25 mm) ID, 1/16th inch (1.6 mm) OD,
3150 mm length, HPLC grade, solid blue, with tubing connectors. ARE-
Applied Research brand. 1 EA.



Additional Info

These Flexible, Coiled PEEK Tubes will expand and contract which will save space in places such as behind an instrument and provide a Pulse Dampener Effect for Parabolic Flow Fluids which will result in a quieter base line and also reduce band broadening due to excess tubing lengths.

This coil is 150 mm wide when resting and 500 mm when fully extended, with a total tubing length of 3150 mm, end to end.

Type 7 FlexiTube, color coded: blue for 0.010 inch (0.25 mm) ID

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$80.89

Printed On 2024/05/15