

## Catalog Number

49330-GC1

**Guard Column Kit, PEEK, 3.0 mm ID, 8 mm active length. ARE-Applied Research brand. 1 EA.**



### Additional Info

This Guard Column Hardware Kit includes one 3.0mm ID Column Tube, two frits and the Guard Column Holder. All components, including the Frits are extruded from ® PEEK. No chemical treatment is done to the tubes and the interior mirror like finish is provided through our very careful extrusion techniques.

Carbon is added to black components.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$126.90**

Printed On 2024/06/20