

Catalog Number

49330-GC2-005

Guard Column Cartridge, PEEK, 3.0 mm ID, 8mm active length. Frits sold separately. ARE-Applied Research brand. 5/PK.



Pack of 5

Additional Info

This Guard Column Cartridge (tube) is 3.0mm ID. Extruded from PEEK, no chemical treatment is done and the interior mirror like finish is provided through our very careful extrusion techniques.

Frits are sold separately.

Price
\$35.79

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$35.79