

## Catalog Number

49330-GC3

**Guard Column Holder Nut, PEEK, for 3.0 mm ID Guard Columns with 8 mm active length. ARE-Applied Research brand. 1 EA.**



### Additional Info

This Replacement, Guard Column Holder Nut is for 3.0mm ID Guard Columns. Extruded from PEEK,

Carbon is added to PEEK to achieve black color.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$31.78**

Printed On 2024/05/05