

## Catalog Number

49330-GC4

**Guard Column Holder Insert**, PEEK, for 3.0 mm ID Guard Columns with 8 mm active length. ARE-Applied Research brand. Bulk Price.



**Bulk Pricing**

### Additional Info

This Replacement, Guard Column Holder Insert is for 3.0mm ID Guard Columns. Extruded from PEEK,

## Additional Info

### **MicroSolv**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)