


Catalog Number

49331-SR2

**Slurry Reservoir Kit, All-PEEK, for packing 4.6 mm ID columns that are longer than 100 mm in length.** Reservoir is 300 mm long with a 7.5 mm ID, complete with all parts and accessories needed to pack bioinert,

## Slurry Reservoir Kit

For 4.6 mm ID columns longer than 100 mm |



1. Slurry Reservoir, Black PEEK 7.5mm ID, 100 mm long
2. PEEK sealing rings / frit spacers
3. *Nothing shown for this kit*
4. End fitting to connect tubing from a pump to reservoir
5. Packing adapter to connect the reservoir to the column

## Additional Info

This kit is for packing All-PEEK columns 4.6 mm ID that are up to 300 mm in length.

The Reservoir Body included in this kit is 300 mm long with 7.5 mm ID and has an internal volume of approximately 13.2 ml. This reservoir is equipped with 5/8-18 external threads on both sides. The adapter has an internal 1/2-20 screw thread on one side and internal 5/8-18 screw threads on the other side. The adapter connects the reservoir to the column being packed.

Kit includes everything you need to pack a bioinert All-PEEK column:

Price

**\$909.45**

## Additional Info

- One Slurry Reservoir - the "body" which is black PEEK.
- One [End Fitting](#) to connect tubing from a pump to the reservoir.
- Three [PEEK sealing rings / frit spacers](#) for the bottom and top connection as shown.
- One [Packing Adapter](#) to connect the reservoir to the column.

How to assemble slurry reservoir and column for use. Click [HERE](#).

**Columns and accessories to be packed are sold separately. Click [HERE](#) for ordering information.**

*This All-PEEK Slurry Reservoir Kit is designed and manufactured for packing only our All-PEEK HPLC columns that are longer than 100 mm with a 4.6 mm ID. This Reservoir and components are matched for maximum performance and are not supported for use with stainless steel or other packing systems. Also, it is recommended to thoroughly clean each part individually and not keep the Packing Reservoir tightly assembled when not in use.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$909.45**