

Catalog Number

49333-EF4

End Fitting, PEEK, for 4.6 mm ID columns, ARE-Applied Research brand. 1 EA.



Additional Info

This End Fitting is made for 4.6mm ID All-PEEK Columns with standard HPLC column 10-32 (tubing inlet) screw threads and 1/2"-20 "column tube" screw threads.

Designed for use only with our All-PEEK HPLC Column Hardware and not recommended for use with stainless steel columns.

Machined from PEEK which ensures a biocompatibility from the beginning of the packing and column use.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$52.96

Printed On 2024/05/02