

## Catalog Number

49333-PS

**Packing Seal, PEEK, 7.5 mm ID, Column Packing Hardware, for use with All-PEEK Column Packing System. ARE-Applied Research brand. 1 EA.**



### Additional Info

This Packing Seal is made for use with the Packing Adapter and 7.5mm ID PEEK Columns. Used for packing All-PEEK HPLC Column Hardware and not recommended for use with stainless steel packing systems.

This Seal is made from PEEK which ensures a biocompatibility from the beginning of the packing and column process.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$24.08**

Printed On 2024/05/02