

## Catalog Number

49407-50

**Tubing, Stainless Steel, Pre-Cut, 0.007 inch (0.18 mm) ID, 1/16 inch (1.60 mm) OD, HPLC grade, black band, 50 cm long. ARE-Applied Research brand. 1 EA.**

**Additional Info**

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut Stainless Steel Tubes are the perfect solution.

Welded and annealed, the ends are cut using Electronic Discharge (*EDC*) Technique to achieve square and burr free zero dead volume connections. They are tested according to EN 10217 7. The Black Band is a color code for .007" ID.

- *Internal and Outer Dimension Tolerances: +/- 0.05mm.*
- *Wall Thickness: 0.706 mm*

**Price**  
**\$29.63**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$29.63**