

## Catalog Number

49410-05

**Tubing, Stainless Steel, Pre-Cut, 0.010 inch (0.254 mm) ID, 1/16 inch (1.59 mm) OD, HPLC Grade, Blue Band, 5cm Long. ARE-Applied Research brand. 1 EA.**



### Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut Stainless Steel Tubes are the perfect solution.

Welded and Annealed, the ends are cut using Electronic Discharge (*EDC*) Technique to achieve square and burr free zero dead volume connections. They are tested according to EN 10204/3.1B. The Blue Band is a color code for .010" ID.

- *Internal and Outer Dimension Tolerances: +/- 0.05mm.*
- Wall Thickness: 0.168 mm

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$20.34**

Printed On 2024/05/04