


Catalog Number

49410-50

Tubing, Stainless Steel, Pre-Cut, 0.010 inch (0.254 mm) ID, 1/16 inch (1.59 mm) OD, HPLC Grade, Blue Band, 50cm Long. ARE-Applied Research brand. 1 EA.



Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut Stainless Steel Tubes are the perfect solution.

Welded and annealed, the ends are cut using Electronic Discharge (*EDC*) Technique to achieve square and burr free zero dead volume connections. They are tested according to EN 10204/3.1B. The Blue Band is a color code for .010" ID.

- *Internal and Outer Dimension Tolerances: +/- 0.05mm.*
- Wall Thickness: 0.168 mm

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$12.43

Printed On 2024/05/08