

## Catalog Number

49420-50

**Tubing, Stainless Steel, Pre-Cut, 0.020 inch (0.508 mm) ID, 1/16 inch (1.59 mm) OD, HPLC Grade, Yellow Band, 50cm Long. ARE-Applied Research brand. 1 EA.**



### Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut Stainless Steel Tubes are the perfect solution.

Welded and annealed, the ends are cut using Electronic Discharge (*EDC*) Technique to achieve square and burr free zero dead volume connections. They are tested according to EN 10204/3.1B. The Yellow Band is a color code for .020" ID.

- *Internal and Outer Dimension Tolerances: +/- 0.05mm.*
- Wall Thickness: 0.541 mm

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$28.98**

Printed On 2024/05/18