

## Catalog Number

49590-05

**Tubing Connector Nuts, high pressure, 1/16 inch, Stainless Steel hex head, medium length with 10-32 threads. Ferrules sold separately. ARE-Applied Research brand. 1 EA.**



### Additional Info

These tubing connector nuts used to connect stainless steel tubing in the high pressures zones of HPLC systems. This 316L stainless steel Fitting has a standard diameter hexagonal head for wrenching the nut.

The medium thread length makes this nut ideal for use on injectors such as the Rheodyne™ and others, where space is limited.

Use with single cone stainless steel ferrules which are designed to work with this fitting, providing a permanent attachment to SS tubing for a high pressure connection.

Click [HERE](#) for ordering information for ferrules.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$7.70**

Printed On 2024/05/02