

Catalog Number

49590-F116-2

Ferrules, stainless steel, two step, single cone, for 1/16 inch tubing.
ARE-Applied Research brand. 10/PK.



Additional Info

These 316L, stainless steel two step ferrules are made to be used with 316L, stainless steel tubing connector, male nuts for connecting 1/16th inch SS tubing to injectors, in high pressures zones of HPLC systems.

The unique shape of this ferrule is designed to be used with Rheodyne™ brand injectors that require ferrules like this.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$113.63

Printed On 2024/04/28