

Catalog Number

58022-N25-LNK

Syringe Filters, 25mm, Nylon, 0.22µm Pore Size. Purple
Polypropylene, AQ Brand. 1500/CS.



Additional Info

Advanced Quality Nylon-66 Syringe Filters, AQ™ brand, for Pharmaceutical Grade HPLC & Dissolution Testing. Manufactured with a unique 0.22µm nylon-66 membrane in a 25mm polypropylene housing, these syringe filters are made for laboratories that require consistent, high performance results with ease of use especially during high volume work. Excellent choice when laboratories cannot afford poor chromatographic results.

These filters are compatible with aqueous or organic HPLC solvents and have been designed to be extremely low adsorbing even for charged molecules and extractable free even in the first 2ml of filtration. The uniquely designed outflow path has 16 channels for minimum back pressure and allows for less hand pressure from the user making for greater compliance and fewer complaints.

Validated for many regulatory requirements, they are designed, tested and QC certified. They comply with US FDA CFR Title 21, Part 210.3(b)(6) for minimum Fiber Release and with USP

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,351.26

Printed On 2024/05/04