

Regenerate Cellulose 4mm Syringe Filters

58022-
R04-K
Part No.

Syringe Filters, 4mm, RC, 0.22µm Pore Size. Grey Polypropylene, AQ Brand. 1000/CS.



Additional Info

Advanced Quality Regenerated Cellulose Syringe Filters, AQ™ brand, for Pharmaceutical Grade HPLC & Dissolution Testing. Manufactured with a unique 0.22µm Regenerated Cellulose Membrane in a 4mm polypropylene housing, these syringe filters are made for laboratories that require consistent, high performance results with ease of use especially during high volume work. Excellent choice when laboratories cannot afford poor chromatographic results.

These filters are used for filtration and sterilization of Protein Solutions, Culture Media and Buffers. The low Protein binding characteristics make these ideal for these types of applications, particularly where the sample size is very small and recovery of protein of interest is sought. Regenerated Cellulose filters provide absolute retention, very low hold up volumes while offering excellent throughputs. Do not use these filters with strong Acids or Bases.

Validated for many regulatory requirements, they are designed, tested and QC certified. They

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,730.25

Printed On 2024/05/01