

## Catalog Number

58022-V25-K

**Syringe Filters, 25mm, PVDF, 0.2 $\mu$ m Pore Size. Dark Purple Polypropylene, AQ Brand. 1000/CS.**  
CUSTOM ITEM



### Additional Info

Advanced Quality Polyvinylidene Fluoride (PVDF) Syringe Filters, AQ™ brand, for Pharmaceutical Grade HPLC & Dissolution Testing. Manufactured with a unique 0.2 $\mu$ m PVDF membrane in a 25mm polypropylene housing, these Syringe Filters are made for Laboratories that require consistent, high performance results with ease of use especially during high volume work. Excellent choice when laboratories cannot afford poor Chromatographic results.

Validated for many Regulatory Requirements, they are designed, tested and QC Certified. They comply with US FDA CFR Title 21, Part 211.92 for minimum Fiber

**Price**

**\$2,698.26**

## Additional Info

Release and with USP <788> test for particulate matter in injectables as well as USP for "Sterile Water for Injection" and therefore do not require any initial flushing or discarding of the first 200µl of filtrate saving time and avoiding investigation of failed methods. The polypropylene plastic housing material complies with the requirements of US FDA CFR Title 21 Part 177.1520 and are heat sealed to prevent leaching. These filters pass a test for oxidizable matter as per the USP.

Color Coded: Dark Purple for 0.2µm PVDF, these filters can be autoclaved for 0.5 hours at 121°C, 0.1MPa or with ethylene oxide. Each Pack is supplied with a Certificate of Analysis.

Click [HERE](#) for 100 in a package quantity or Click [HERE](#) for 50 in a package quantity.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,698.26**