

## Catalog Number

## 58045-NDX-100

**Syringe Filters, 25mm, Nylon, w/ Multi-Depth Pre Filters,**  
**0.45µm** Final Pore Size. Black with White Polypropylene, AQ Brand.  
100/PK.



### Additional Info

Advanced Quality Nylon-66 Syringe Filters with 3 Binder Free, Borosilicate, Micro-glass-fiber "Pre-Filters", AQ™ brand, for Pharmaceutical Grade HPLC & Dissolution Testing where samples need multiple layers of pre-filters. Manufactured with a unique 0.45µm nylon-66 membrane in a 25mm polypropylene housing, these syringe filters are made for laboratories that require consistent, high performance results with ease of use especially during high volume work. Excellent choice when laboratories cannot afford poor results.

These filters are compatible with aqueous or organic HPLC solvents and have been designed for high through put and to be extractable free even in the first 2ml of filtration. The uniquely designed outflow path has 16 channels for minimum back pressure and allows for less hand pressure from the user making for greater compliance and fewer complaints.

Validated for many regulatory requirements, they are designed, tested and QC certified. They comply with US FDA CFR Title 21, Part 210.3(b)(6) for minimum Fiber Release and with USP

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$415.83**

Printed On 2024/04/25