

## Regenerated Cellulose 4mm Syringe Filter

**58045-  
R04-K**  
Part No.

**Syringe Filters, 4mm, RC, 0.45um** Pore Size. Brown Polypropylene, AQ Brand. 1000/CS.



### Additional Info

Advanced Quality Regenerated Cellulose Syringe Filters, AQ™ brand, for Pharmaceutical Grade HPLC & Dissolution Testing. Manufactured with a unique 0.45µm Regenerated Cellulose Membrane in a 4mm polypropylene housing, these syringe filters are made for laboratories that require consistent, high performance results with ease of use especially during high volume work. Excellent choice when laboratories cannot afford poor chromatographic results.

These filters are used for filtration and sterilization of Protein Solutions, Culture Media and Buffers. The low Protein binding characteristics make these ideal for these types of applications, particularly where the sample size is very small and recovery of protein of interest is sought. Regenerated Cellulose filters provide absolute retention, very low hold up volumes while offering excellent throughputs. Do not use these filters with strong Acids or Bases.

Validated for many regulatory requirements, they are designed, tested and QC certified. They

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,730.25**

Printed On 2024/05/02