

## Catalog Number

58210-22

**Needles, Blunt Tip 90 degrees**, 20 gauge x 1.5 inches, Non Sterile for Laboratory Use Only. BASIK Brand. 1000/PK



### Additional Info

These non-sterile needles can be used with the Luer Lock or Luer Slip, MicroSolv Chrom Syringes in the laboratory. The needle is stainless steel and the hub is plastic material. The syringes come with a protective shield.

The 90° blunt tip is good for syphoning from vials, bottles, biological systems or any small pool of liquid. These syringes can be connected to plastic tubing to act a small, low pressure "pump".

*Not recommended for piercing autosampler vial septa.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$215.73**

Printed On 2024/04/28