

Catalog Number

58800-47-KT-2L

Mobile Phase Filter Degasser Assembly, Glass, 300ml Reservoir, 2L Flask, Complete. BASIK Brand. 1 EA.



Additional Info

This borosilicate glass, Mobile Phase Filtration Assembly is great for filtering and degassing HPLC solvents for use as mobile phases.

The unit's design of the fritted glass support for the membrane also prevents contamination of your vacuum line or source which is separate from this unit. The funnels are graduated and will hold up to 300mL of mobile phase. The reservoir flasks holds 2,000mL and the narrow mouth opening minimizes out-gassing. Use any disposable, 47mm membrane disc that is compatible with your mobile phase. The Filter Degasser Assembly is sold as a complete kit including the blue metal clamp. Replacements parts are available upon request. *Membranes are sold separately.*

Use caution when selecting a vacuum pump that is appropriate for glass filter devices.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$488.46

Printed On 2024/05/04