

Catalog Number

58801-02

Bottle, Mobile Phase Reservoir, Safety, Coated Glass, Clear, 1L,
Standard Taper Top, 29/32 Neck. MicroSolv Brand. 1 EA.



Additional Info

This Plastic Coated, Borosilicate Glass Bottle is perfect for HPLC Mobile Phase Reservoirs with its small Opening and Large Capacity. Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle requires a S/T (Standard Taper) 29/32 joint for the Mouth and will produce a tight seal. These Bottles are often used on older HPLC Instruments such as the Agilent 1090 and with Filter Degasser similar to [58800-00](#) (sold separately) or other Degassers.

The outside only of the Bottle is coated in a transparent, thick epoxy plastic and if this bottle is dropped and broken, the plastic will prevent a spill incident and prevent workers from getting cut with glass or splashed with solvent.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$223.98

Printed On 2024/04/28