

Catalog Number

58802-05

Bottle, Mobile Phase Reservoir, Safety, Coated Glass, Clear, 5L,
Screw Top, GL45 Threads. MICROSOLV Brand. 1 EA.



Additional Info

This Clear, Plastic Coated, Borosilicate Glass Bottle has a 5L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir. The small opening, Large Capacity and Enamel Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MICROSOLV Filter Degassers](#), other Degassers or disposable [Filter Degassers](#). A Solid Screw Cap without a liner and a polypropylene Drip Ring are included.

Price
\$620.49

Additional Info

Larger bottles like this one are often used for making mobile phase at the start of the day and throughout the day the bottles on the instrument can be topped off with the same mobile phase minimizing issues from mobile phase preparation. This is especially useful for overnight runs before the end of the workday.

This Bottle is often used with the [Mobile Phase Delivery System](#). The outside only of the Bottle is coated in a transparent, thick epoxy plastic and if this bottle is dropped and broken, the plastic will prevent a spill incident and prevent workers from getting cut with glass or splashed with solvent.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$620.49