

## Catalog Number

58802-20-R

Cap, Screw Top GL45 Threads, Red, Open Hole. MicroSolv Brand. 1 EA.



### Additional Info

This cap is made of Polypropylene and has GL45 threads making it a good choice for Mobile Phase Reservoir Bottles.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$15.29**

Printed On 2024/05/18