

Catalog Number

58802-20-R

Cap, Screw Top GL45 Threads, Red, Open Hole. MicroSolv Brand. Bulk Price.



Bulk Pricing

Additional Info

This cap is made of Polypropylene and has GL45 threads making it a good choice for Mobile Phase Reservoir Bottles.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/07/27