

## Catalog Number

58807-050-CS

**Bottles, HPLC Mobile Phase Storage**, Glass, Clear, 50ml, Screw Top, GL32 Threads. MicroSolv Brand. 10/CS.



### Additional Info

These Clear, Borosilicate Glass Bottles have a 50ml "Total Fill Capacity" and are perfect as HPLC Mobile Phase Storage Bottles. The small openings and Enamel Graduations and Marking Patches makes them a convenient way to store concentrated Mobile Phase or its Ingredients.

Blue, Solid Screw Caps without liners and polypropylene Drip Rings are included. These Bottle are *Not "Safety Coated"*.

Click [HERE](#) for smaller package quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$226.75**

Printed On 2024/05/03