

Catalog Number

58807-100-A

Bottle, HPLC Mobile Phase Storage, Glass, Amber, 100ml, Screw Top, GL45 Threads. MicroSolv Brand. 1 EA.



Additional Info

This Amber, Borosilicate Glass Bottle has a 100ml "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Storage Bottle. The small opening and Enamel Graduations and Marking Patch makes this a convenient way to store concentrated Mobile Phase or its Ingredients.

A Blue, Solid Screw Cap without a liner and a polypropylene Drip Ring is included. This Bottle is *Not "Safety Coated"*.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$93.80

Printed On 2024/05/05