

## Catalog Number

58807-250

**Bottle, HPLC Mobile Phase Storage**, Glass, Clear, 250ml, Screw Top, GL45 Threads. MicroSolv Brand. 1 EA.



### Additional Info

This Clear, Borosilicate Glass Bottle has a 250ml "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Storage Bottle. The small opening and Enamel Graduations and Marking Patch makes this a convenient way to store concentrated Mobile Phase or its Ingredients.

A Blue, Solid Screw Cap without a liner and a polypropylene Drip Ring is included. This Bottle is *Not "Safety Coated"*.

Click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$25.88**

Printed On 2024/04/29