

Catalog Number

58807-500

Bottle, HPLC Mobile Phase Storage, Glass, Clear, 500ml, Screw Top, GL45 Threads. MICROSOLV Brand. 1 EA.



Additional Info

This Clear, Borosilicate Glass Bottle with a 500 mL "Total Fill Capacity" is an essential tool for any laboratory. Its durable Hybex™ borosilicate glass construction and precise white enamel graduations make it perfect for use as an HPLC mobile phase storage bottle.

- Material: Made from high-quality Hybex™ borosilicate glass (3.3), known for its superior chemical resistance and durability.
- Capacity: The bottle has a total fill capacity of 500 milliliters, making it ideal for

Price
\$34.24



Additional Info

use as an HPLC mobile phase storage.

- Design: Features a small opening that helps minimize evaporation and contamination. The white enamel graduations on the outside allow for easy monitoring of solvent levels, ensuring you don't run dry.
- Temperature Range: Can withstand temperatures from -70°C to 500°C, making it suitable for use in freezers, microwaves, autoclaves, and dry heat sterilization chambers.
- Cap: Comes with a standard blue polypropylene GL45 sealing cap, which is autoclavable and can withstand temperatures up to 140°C and a polypropylene Drip Ring.

This Bottle is *Not "Safety Coated"*.

Click [HERE](#) for case quantities.

MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$34.24