

Catalog Number

58807-Z32

Pouring Ring, Replacement, Natural, for HPLC Mobile Phase Reservoir Bottles. GL32 Threads. MicroSolv Brand. 10/PK.



Additional Info

Using this Natural color Polypropylene Pouring or "Drip Ring" can make decanting liquids from any bottle with a GL32 Thread easier and prevent that last drop from getting away from you.

This ring will fit deeper into the mouth of the Bottle, making is more secure and less likely to fall out.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$22.50

Printed On 2024/07/20