

## Catalog Number

## 58819-05R

**Bottle, Mobile Phase Reservoir, Safety, Coated Glass, Clear, 2L, Screw Top, GL45 Threads, Open Red Cap. MicroSolv Brand. 1 EA.**



### Additional Info

This Clear, Plastic Coated, Borosilicate Glass Bottle has a 2L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir. The small opening, two-liter Capacity and White Enamel Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MicroSolv Filter Degassers](#), other Degassers, disposable [Filter Degassers](#). A Red, Open Top Screw Cap and a Red polypropylene Pouring Ring are included.

*This Bottle is often used with the [Mobile Phase Delivery System](#). The outside only of the Bottle is coated in a transparent, thick epoxy plastic and if this bottle is dropped and broken, the plastic will prevent a spill incident and prevent workers from getting cut with glass or splashed with solvent.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$255.53**

Printed On 2024/07/20